

NEW KU-3™ Viscometer

for Paints, Coatings, and Inks

ASTM D562 Compatible
(industry specification)

New Magnetic Spindle Coupling

New User Interface
with touch control & auto-start

New Single Piece Can Adapter
for pint and 1/2 pint cans
Accommodates quart cans

Easy to use
no weights, simplifies an
established test procedure

LED Display Info:

- Krebs Units
- Gram Units
(Weight)
- Centipoise*

Select Krebs or Grams or Centipoise

Lock-In Test Results
with Hold button

Accuracy: $\pm 1.0\%$ of range

Repeatability: $\pm 0.5\%$

Standard Krebs Spindle

Measurement range:

40 to 141 KU, 32 to 1099 gm,
and 27 to 5274 cP*



What's Included?

Instrument
Krebs-type Spindle (p48)
Adapter for Quart, Pint and
Half-Pint Cans

Optional Accessories

Paste Spindle (p48)

Applications

Paints
Coatings
Adhesives
Inks
Pastes

*Centipoise values based on the conversion from Krebs
Units as defined in the ASTM D562.