








MEASUREMENT PARAMETERS

MODEL	 WEIGHT	 THICKNESS	 DIAMETER/ LENGTH	 HARDNESS	 WIDTH
P2	●	●	—	—	—
P3	—	●	●	●	—
P4	●	●	●	●	—
P5	●	●	●	●	●

MODEL	W x D x H (mm)	Weight (kg)
P2 - P5	320 x 320 x 185	< 15

Tablet positioning that's one of a kind with the Vibralign and Roto systems

The P5-series sets new standards in the laboratory testing technology. The P5 version tests all five parameters automatically. The Vibralign and Roto systems make it easy to position tablets precisely.



Vibralign system

The tablets are precisely rotated 90° and positioned for width measurement using the Vibralign system incorporated in the hardness station.



Roto system

(manual, adjustable grooved plate)

In the P4 model, the Roto system (manual, adjustable grooved plate) ensures easy tablet positioning. You can choose from three different settings: a flat, light or deep groove.

ADVANTAGES

- + Sensor-aided navigation
- + Integrated 360° LED status display
- + Automatic tablet positioning
- + Clean design

OPTIONS

- Load cells: 50 N, 500 N (standard), 1000 N
- Covered model with integrated breakage chamber
- Vibralign system
- Roto system (manual, adjustable grooved plate) for simple positioning of test specimens in the P3 and P4 models



The P-series always provides the best possible solution for your requirements.



Intuitive navigation

- + Insert tablets and start testing
- + The new touchscreen user interface of the P-series is simple and self-explanatory



Integrated 360° LED status display

- + The clear LED display provides information on the unit's status at all times
- + The status of the unit can easily be seen even from a distance



Automatic tablet positioning

- + Length, width, hardness, ...
- + The P-series accurately places even oblong tablets in every test position



Clean design

- + Always hygienic and clean
- + With its rounded corners and edges, the P-series is easy to clean